


Description

Product Name	Norethynodrel
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor agonist
Calculated MW	298.42
Formulation	C ₂₀ H ₂₆ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@]2(C#C)O)CCC4=C3CCC(=O)C4Pathways:Endocrinology/HormonesReceptor:ERBolin
g pt:441.6CMelting pt:169-170CSolubility:Water: 8.95 mg/L .Appearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Edgren RA. Steroids. 1994 Jan;59(1):58-9.

Note: This product is for in vitro research use only