

Nylidrin hydrochloride

Catalog No: #S0173


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nylidrin hydrochloride
Brief Description	Inhibitors
Purification	99.50%
Target Name	Adrenergic Receptor agonist
Calculated MW	335.87
Formulation	C19H26ClNO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: Cl.CC(CCC1=CC=CC=C1)NC(C)C(O)C1=CC=C(O)C=C1
 Pathways: GPCR/G Protein Receptor; β -adrenergic receptor
 Boiling pt: Melting pt: 198-199°C
 Solubility: DMSO: 38mg/mL
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Garetz FK, et al. J Am Geriatr Soc. 1979 May;27(5):235-6.

Note: This product is for in vitro research use only