

(S)-Clopidogrel

Catalog No: #S0182


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	(S)-Clopidogrel
Brief Description	Inhibitors
Purification	99.50%
Target Name	P2 Receptor antagonist
Calculated MW	321.82
Formulation	C16H16ClNO2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: COC(=O)[C@H](C1=CC=CC=C1Cl)N2CCC3=C(C2)C=CS3 Pathways: Neuroscience Receptor: P2Y12 Boling pt: 423.7°C Melting

pt: Solubility: Water: 50.78 mg/L Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Dorsam RT, et al. Hematology. 2003 Dec;8(6):359-65.

Note: This product is for in vitro research use only