

() Clopidogrel hydrogensulfate

Catalog No: #S0182L2


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	() Clopidogrel hydrogensulfate
Brief Description	Inhibitors
Purification	98.00%
Target Name	P2 Receptor antagonist
Calculated MW	419.9
Formulation	C16H16ClNO2S H2SO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: COC(=O)[C@H](C1=CC=CC=C1Cl)N2CCC3=C(C2)C=CS3.OS(=O)(=O)O Pathways: Neuroscience Receptor: P2Y12

Bolting pt: 423.7°C Melting pt: 184°C Solubility: DMSO: 10 mM Appearance: Colorless or Faint Yellow Liquid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Dorsam RT, et al. Hematology. 2003 Dec;8(6):359-65.

Note: This product is for in vitro research use only