

## Pramiracetam

Catalog No: #S0183


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pramiracetam
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	269.38
Formulation	C <sub>14</sub> H <sub>27</sub> N <sub>3</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
SMILES: CC(C)N(CCNC(=O)CN1CCCC1=O)C(C)C  
Pathways: Others  
Receptor: Others  
Boiling pt: 461°C  
Melting pt: 100°C  
Solubility: DMSO: 54 mg/mL (200.46 mM)  
Appearance: White Crystal  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
Stock solution can be stored below -20 °C for several month.

## References

1. Corasaniti MT, et al. *Funct Neurol*. 1995 May-Jun;10(3):151-5.

Note: This product is for in vitro research use only