

## Secnidazole

Catalog No: #S0195


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Secnidazole
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	185.18
Formulation	C7H11N3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: CC(O)CN1C(C)=NC=C1[N+](=O)[O-] Pathways: Others Receptor: Others Boiling pt: Melting pt: 76°C Solubility: Ethanol: 37 mg/mL (199.8 mM); Water: 37 mg/mL (199.8 mM); DMSO: 37 mg/mL (199.8 mM) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. De Backer, et al. Clinical Microbiology and Infection 16 (5): 470C472.

Note: This product is for in vitro research use only