

Sulbentine

Catalog No: #S0202


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sulbentine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	314.47
Formulation	C17H18N2S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:S=C1SCN(CC2=CC=CC=C2)CN1CC1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:451.5CMelting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Thoma K, et alPharmazie. 1997 May;52(5):362-73.

Note: This product is for in vitro research use only