

Protriptyline hydrochloride

Catalog No: #S0204


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protriptyline hydrochloride
Brief Description	Inhibitors
Purification	98.54%
Target Name	5-HT Receptor inhibitor
Calculated MW	299.84
Formulation	C19H21NHCl
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:Cl.CNCCCC1C2=CC=CC=C2C=CC2=CC=CC=C12Pathways:GPCR/G Protein; NeuroscienceReceptor:Sert (Sodium-dependent)Boiling pt:407.7CMelting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. McDougale CJ, et al. Mol Psychiatry. 1998 May;3(3):270-3.

Note: This product is for in vitro research use only