

Prednisolone phosphate sodium

Catalog No: #S0208


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prednisolone phosphate sodium
Brief Description	Inhibitors
Purification	99.50%
Target Name	GR inhibitor
Calculated MW	484.39
Formulation	C ₂₁ H ₂₇ Na ₂ O ₈ P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12C[C@@H]([C@H]3[C@H]([C@@H]1CC[C@@]2(C(=O)COP(=O)([O-])[O-])O)CCC4=CC(=O)C=C[C@]34C)O.[Na+].[Na+]Pa

thways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:100CMelting pt:Solubility:Soluble in WaterAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ikonomidis I, et al. Heart. 2009 Sep;95(18):1502-7.

Note: This product is for in vitro research use only