

Azelastine

Catalog No: #S0211


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Azelastine
Brief Description	Inhibitors
Purification	99.50%
Target Name	Histamine Receptor antagonist
Calculated MW	381.9
Formulation	C22H24ClN3O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CN1CCCC(CC1)N1N=C(CC2=CC=C(Cl)C=C2)C2=CC=CC=C2C1=O
 Pathways: GPCR/G Protein;
 Immunology/Inflammation; Neuroscience
 Receptor: H1 receptor
 Boiling pt: 533.9°C
 Melting pt: Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Casale TB. J Allergy Clin Immunol. 1989 Apr;83(4):771-6.

Note: This product is for in vitro research use only