

CLIC4 Polyclonal Antibody

Catalog No: #40758

Package Size: #40758-1 50ul #40758-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	CLIC4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	CLIC4 Polyclonal Antibody detects endogenous levels of CLIC4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human CLIC4.
Target Name	CLIC4
Other Names	CLIC4; Chloride intracellular channel protein 4; Intracellular chloride ion channel protein p64H1
Accession No.	Swiss-Prot: Q9Y696NCBI Gene ID: 25932
Uniprot	Q9Y696
GeneID	25932;
SDS-PAGE MW	29kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of HT29 cells using CLIC4 Polyclonal Antibody



Western Blot analysis of L929 RAT-MUScle rat-kidney cells using CLIC4 Polyclonal Antibody

Note: This product is for in vitro research use only