

Seratrodast

Catalog No: #S0212


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Seratrodast
Brief Description	Inhibitors
Purification	98.93%
Target Name	COX inhibitor; PPAR antagonist
Calculated MW	354.44
Formulation	C22H26O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory system
 SMILES:CC1=C(C)C(=O)C(=C(C)C1=O)C1=CC(CCCCCC(O)=O)=CC=C1
 Pathways:Immunology/Inflammation; Metabolism; Neuroscience
 Receptor:COX; TXA2R
 Boling pt:522.3C?Melting pt:118-120C
 Solubility:Soluble in Water
 Appearance:Yellow Powder
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Endo S, Akiyama K. Nihon Rinsho. 1996 Nov;54(11):3045-8.

Note: This product is for in vitro research use only