

## Pioglitazone

Catalog No: #S0214


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Pioglitazone
Brief Description	Inhibitors
Purification	90.00%
Target Name	PPAR agonist
Calculated MW	356.44
Formulation	C19H20N2O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism

systemSMILES:Cl.CCC1=CN=C(CCOC2=CC=C(CC3SC(=O)NC3=O)C=C2)C=C1Pathways:MetabolismReceptor:PPARBoling pt:575.4CMelting pt:183-184CSolubility:DMSO: 15 mg/mL (42.08 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Jaakkola T, et al. Basic Clin Pharmacol Toxicol. 2006 Jul;99(1):44-51.

Note: This product is for in vitro research use only