

# Losartan potassium


Catalog No: #S0215

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Losartan potassium
Brief Description	Inhibitors
Purification	99.50%
Target Name	RAAS antagonist
Calculated MW	461.01
Formulation	C22H22ClKN6O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular

System SMILES: [K].CCCCC1=NC(Cl)=C(CO)N1CC1=CC=C(C=C1)C1=CC=CC=C1C1=NN=N[N-]1 Pathways: Endocrinology/Hormones Receptor: AT1 receptor  
 Boiling pt: 134°C Melting pt: 76.9°C Solubility: Ethanol: 92 mg/mL (199.12 mM); Water: 92 mg/mL (199.12 mM); DMSO: 92 mg/mL (199.12 mM)  
 Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Burnier M, et al. Circulation, 2001, 103, 904-912.

Note: This product is for in vitro research use only