

Losartan

Catalog No: #S0215L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Losartan
Brief Description	Inhibitors
Purification	99.48%
Target Name	RAAS antagonist
Calculated MW	422.91
Formulation	C22H23ClN6O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: CCCCC1=NC(Cl)=C(CO)N1CC=CC=C(C=C1)C1=CC=CC=C1C1=NNN=N1 Pathways: Endocrinology/Hormones Receptor: AT1

receptor Boiling pt: 686.5°C Melting pt: 183.5-184.5°C Solubility: Water: 0.82 mg/L Appearance: Faint Yellow Solid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Burnier M, et al. Circulation, 2001, 103, 904-912.

Note: This product is for in vitro research use only