

## Eletriptan HBr

Catalog No: #S0216


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Eletriptan HBr
Brief Description	Inhibitors
Purification	99.50%
Target Name	5-HT Receptor agonist
Calculated MW	463.43
Formulation	C22H27BrN2O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: CN1CCC[C@@H]1CC2=CNC3=C2C=C(C=C3)CCS(=O)(=O)C4=CC=CC=C4.Br  
 Pathways: GPCR/G Protein;  
 Neuroscience Receptor: 5-HT1B; 5-HT1DB  
 Boiling pt: 633.9C  
 Melting pt: 169-171C  
 Solubility: DMSO: 93 mg/mL (200.67 mM)  
 Appearance: Off-white  
 Color Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Napier C, et al. Eur J Pharmacol. 1999 Mar 5;368(2-3):259-68.

Note: This product is for in vitro research use only