

## Foscarnet sodium

Catalog No: #S0220


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Foscarnet sodium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Reverse Transcriptase inhibitor
Calculated MW	191.95
Formulation	CNa3O5P
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: C(=O)([O-])P(=O)([O-])[O-].[Na+].[Na+].[Na+]  
 Pathways: Microbiology&virology  
 Receptor: Reverse Transcriptase;  
 Reverse Transcriptase  
 Boiling pt: 270-272°C  
 Melting pt: 62-65°C  
 Solubility: Soluble in DMSO  
 Appearance: White Crystal  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- Oberg B. Pharmacol Ther, 1989, 40(2), 213-285.

Note: This product is for in vitro research use only