

Thioproline

Catalog No: #S0221


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thioproline
Brief Description	Inhibitors
Purification	99.00%
Target Name	Annexin A antagonist
Calculated MW	133.17
Formulation	C ₄ H ₇ NO ₂ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: OC(=O)C1CSCN1
Pathways: GPCR/G Protein Receptor: Annexin A5
Boiling pt: 350.3°C
Melting pt: 199-201°C
Solubility: DMSO: 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Barba-Espin G, et al. Plant Physiol Biochem. 2012 Oct;59:30-6.

Note: This product is for in vitro research use only