

Tetroquinone

Catalog No: #S0223


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tetroquinone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	172.09
Formulation	C6H4O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:OC1=C(O)C(=O)C(O)=C(O)C1=OPathways:OthersReceptor:OthersBoiling pt:148CMelting pt:300CSolubility:DMSO: 10 mMAppearance:Black SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lou K, et al. J Am Chem Soc. 2010 Dec 15; 132(49):17635-41.

Note: This product is for in vitro research use only