

Pizotifen Malate

Catalog No: #S0225


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pizotifen Malate
Brief Description	Inhibitors
Purification	99.50%
Target Name	5-HT Receptor agonist
Calculated MW	429.53
Formulation	C19H21NSC4H6O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: OC(=O)C=CC(O)=O.CN1CCC(CC1)=C1C2=C(CCC3=CC=CC=C13)SC=C2
 Pathways: GPCR/G Protein;
 Neuroscience Receptor: 5-HT
 Boiling pt: 436.7°C
 Melting pt: Solubility: DMSO: 21 mg/mL (48.89 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Margarit MV, et al. Drug Dev Ind Pharm. 2001 Jul;27(6):517-22.

Note: This product is for in vitro research use only