

Betaxolol hydrochloride

Catalog No: #S0226


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Betaxolol hydrochloride
Brief Description	Inhibitors
Purification	99.56%
Target Name	Adrenergic Receptor antagonist
Calculated MW	343.89
Formulation	C18H30ClNO3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System SMILES: Cl.CC(C)NCC(O)COC1=CC=C(CCOCC2CC2)C=C1 Pathways: GPCR/G

Protein Receptor: β 1-adrenergic receptor Boiling pt: 448°C Melting pt: Solubility: Ethanol: 69 mg/mL (200.64 mM); Water: 69 mg/mL (200.64 mM); DMSO: 69 mg/mL (200.64 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Tang LQ, et al. Curr Eye Res, 1998, 17(1), 24-30.

Note: This product is for in vitro research use only