

Rivastigmine

Catalog No: #S0229


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rivastigmine
Brief Description	Inhibitors
Purification	99.85%
Target Name	AChR inhibitor
Calculated MW	250.34
Formulation	C14H22N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CCN(C)C(=O)OC1=CC=CC(=C1)[C@H](C)N(C)C
 Pathways: Neuroscience
 Receptor: AChE
 Boiling pt: 316.2°C
 Melting pt: Solubility: Soluble in DMSO
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Kennedy JS, et al. J Clin Psychopharmacol. 1999 Dec;19(6):513-21.

Note: This product is for in vitro research use only