

CLN5 Polyclonal Antibody

Catalog No: #40760

Package Size: #40760-1 50ul #40760-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CLN5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	CLN5 Polyclonal Antibody detects endogenous levels of CLN5 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CLN5.
Target Name	CLN5
Other Names	CLN5; Ceroid-lipofuscinosis neuronal protein 5; Protein CLN5
Accession No.	Swiss-Prot: O75503NCBI Gene ID: 1203
Uniprot	O75503
GeneID	1203;
SDS-PAGE MW	48kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T cells using CLN5 Polyclonal Antibody

Note: This product is for in vitro research use only