

Ethamivan


Catalog No: #S0233

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ethamivan
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	223.27
Formulation	C ₁₂ H ₁₇ NO ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:CCN(CC)C(=O)C1=CC(OC)=C(O)C=C1Pathways:OthersReceptor:OthersBoiling pt:401.4CMelting pt:6-97CSolubility:Water: 12.4 mg/mLAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Hirsh K, Wang SC. J Pharmacol Exp Ther. 1975 May; 193(2):657-63.

Note: This product is for in vitro research use only