

Modaline sulfate

Catalog No: #S0235


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Modaline sulfate
Brief Description	Inhibitors
Purification	99.65%
Target Name	MAO inhibitor
Calculated MW	275.33
Formulation	C10H17N3O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: OS(O)(=O)=O.CC1=NC=CN=C1N1CCCCC1
 Pathways: Metabolism
 Receptor: MAO
 Boiling pt: 294.2°C
 Melting pt: Solubility: Water: 10 mM
 Appearance: Black Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Tangri KK, et al. Jpn J Pharmacol. 1972 Oct;22(5):599-604.

Note: This product is for in vitro research use only