

## Closantel

Catalog No: #S0237


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Closantel
Brief Description	Inhibitors
Purification	97.00%
Target Name	Others
Calculated MW	663.07
Formulation	C22H14Cl2I2N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CC1=CC(C(C#N)C2=CC=C(CI)C=C2)=C(CI)C=C1NC(=O)C1=CC(I)=CC(I)=C1O Pathways: Others Receptor: Antiparasitic Boli  
 ng pt: Melting pt: Solubility: DMSO: 100 mg/mL (150.81 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Hlasta DJ, et al. Bioorg Med Chem Lett, 1998, 8(14), 1923-1928.

Note: This product is for in vitro research use only