

# 2-(Acetyloxy)-3-Methylbenzoic Acid


Catalog No: #S0238

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	2-(Acetyloxy)-3-Methylbenzoic Acid
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others inhibitor
Calculated MW	194.18
Formulation	C10H10O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CC(=O)OC1=C(C=CC=C1C)C(=O)OPathways:OthersReceptor:OthersBoiling pt:324.5CMelting pt:115CSolubility:DMSO: 10 mMAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

N/A

Note: This product is for in vitro research use only