

Vanitiolide


Catalog No: #S0243

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vanitiolide
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	253.32
Formulation	C ₁₂ H ₁₅ NO ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: COC1=CC(=C(N2CCOCC2)S)/C=CC1=O
 Pathways: Others
 Receptor: Others
 Boiling pt: 392.9°C
 Melting pt: 100°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. QHTS Assay For Inhibitors Of Aldehyde Dehydrogenase 1 (ALDH1A1)

Note: This product is for in vitro research use only