

Tilorone dihydrochloride

Catalog No: #S0244


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|-------------------------------------------------------------------------------|
| Product Name | Tilorone dihydrochloride |
| Brief Description | Inhibitors |
| Purification | 97.00% |
| Target Name | Others |
| Calculated MW | 483.47 |
| Formulation | C ₂₅ H ₃₆ Cl ₂ N ₂ O ₃ |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research

Area: Immunology
 SMILES: [H+].[H+].[Cl-].[Cl-].CCN(CC)CCOC1=CC2=C(C=C1)C1=C(C=C(OCCN(CC)CC)C=C1)C2=O
 Pathways: Others
 Receptor: interferon inducer
 Boiling pt: 552.3°C
 Melting pt: 230-233°C
 Solubility: Water: 10 mM
 Appearance: Orange Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

- Stringfellow DA, Glasgow LA. Antimicrob Agents Chem Others. 1972 Aug;2(2):73-8.

Note: This product is for in vitro research use only