

Mycophenolate mofetil

Catalog No: #S0245


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mycophenolate mofetil
Brief Description	Inhibitors
Purification	99.00%
Target Name	Dehydrogenase inhibitor
Calculated MW	433.5
Formulation	C23H31NO7
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism

systemSMILES:CC1=C(C(=C(C2=C1COC2=O)O)C)/C=C(\C)/CCC(=O)OCCN3CCOCC3)OCP

Pathways:MetabolismReceptor:IMPDHBoiling pt:637.6CMelting pt:93-94CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Nakanishi T, et al. Int Immunopharmacol. 2010, 10(1), 91-97.

Note: This product is for in vitro research use only