

## Acecainide


Catalog No: #S0249

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Acecainide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Sodium Channel
Calculated MW	277.36
Formulation	C15H23N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular System  
SMILES: CCN(CC)CCNC(=O)C1=CC=C(NC(C)=O)C=C1  
Pathways: Membrane transporter/Ion channel  
Receptor: Sodium Channel  
Boiling pt: 500°C  
Melting pt: 138-140°C  
Solubility: Soluble in DMSO  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Weiss R, et al. Circulation. 2002 Feb 12;105(6):707-13.

Note: This product is for in vitro research use only