

Mepiroxol

Catalog No: #S0257


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mepiroxol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	125.13
Formulation	C6H7NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OCC1=C[N+]([O-])=CC=C1
Pathways: Others
Receptor: Others
Boiling pt: 391.3°C
Melting pt: 87-89°C
Solubility: Soluble in DMSO
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Robert A, et al. Dig Dis Sci. 1987 Sep;32(9):997-1003.

Note: This product is for in vitro research use only