

Carsalame

Catalog No: #S0258


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Carsalame
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	163.13
Formulation	C8H5NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:O=C1NC(=O)C2=CC=CC=C2O1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:228-232CSolubility:Soluble in DMSOAppearance:White CrystalRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kamal A. Biochem Pharmacol. 1990 Oct 1;40(7):1669-71.

Note: This product is for in vitro research use only