

## Merbromin

Catalog No: #S0265


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Merbromin
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	750.65
Formulation	C <sub>20</sub> H <sub>8</sub> Br <sub>2</sub> HgNa <sub>2</sub> O <sub>6</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: O.[Na+].OC(=O)C1=CC=CC=C1C1=C2C=C(Br)C(=O)C=C2OC2=C([Hg])C(O)=C(Br)C=C12 Pathways: Others Receptor: Others Boli  
 ng pt: Melting pt: 300°C Solubility: DMSO: 10 mM Appearance: Black Solid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Banerjee D, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2008 Oct;70(5):1109-13.

Note: This product is for in vitro research use only