

Dichlorisone Acetate

Catalog No: #S0269


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dichlorisone Acetate
Brief Description	Inhibitors
Purification	99.00%
Target Name	AChR inhibitor
Calculated MW	455.4
Formulation	C ₂₃ H ₂₈ Cl ₂ O ₅
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: CC(=O)OCC(=O)[C@]1(CC[C@@H]2[C@@]1(C[C@@H]([C@]3([C@H]2CCC4=CC(=O)C=C[C@@]43)Cl)Cl)C)O
 Pathways: Neuroscience
 Receptor: AChE
 Boiling pt: 596.2°C
 Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. JEANDIDIER P. Ann Med Nancy. 1964 Aug-Sep;42:1062-5.

Note: This product is for in vitro research use only