

Hydrocortisone butyrate

Catalog No: #S0270


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hydrocortisone butyrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Annexin A; GR
Calculated MW	432.55
Formulation	C25H36O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:CCCC(=O)O[C@@]1(CC[C@@H]2[C@@]1(CC[C@@H]([C@@H]3[C@H]2CCC4=CC(=O)CC[C@]34C)O)C(=O)COPathways:GPCR/
G Protein; MetabolismReceptor:Annexin A1; Glucocorticoid ReceptorBoiling pt:585.6CMelting pt:Solubility:Soluble in WaterAppearance:White
PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be
stored below -20 °C for several month.

References

1. Grossman R, et al. Ann N Y Acad Sci. 2006 Jul;1071:410-21.

Note: This product is for in vitro research use only