

Iodipamide


Catalog No: #S0271

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Iodipamide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	1139.76
Formulation	C ₂₀ H ₁₄ I ₆ N ₂ O ₆
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Endocrine

system SMILES: OC(=O)C1=C(I)C(NC(=O)CCCCC(=O)NC2=C(I)C=C(I)C(C(O)=O)=C2I)=C(I)C=C1I Pathways: Others Receptor: Others Boiling pt: 908.6°C Melting pt: 306-308°C Solubility: Water: 460 mg/L (at 20°C) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Fidler JL, et al. AJR Am J Roentgenol. 2013 Aug; 201(2):W271-82.

Note: This product is for in vitro research use only