

Mefexamide hydrochloride

Catalog No: #S0273


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mefexamide hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	316.82
Formulation	C ₁₅ H ₂₅ ClN ₂ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: [H+].[Cl-].CCN(CC)CCNC(=O)COC1=CC=C(OC)C=C1
 Pathways: Others
 Receptor: Others
 Boling pt: 460.4°C
 Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

1. Vinogradov VM, Grechko AT. Farmakol Toksikol. 1982 Sep-Oct;45(5):108-11.

Note: This product is for in vitro research use only