

## Meticrane

Catalog No: #S0275


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Meticrane
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	275.34
Formulation	C10H13NO4S2
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Other  
 SMILES:CC1=C(C=C2C(CCCS2(=O)=O)=C1)S(N)(=O)=O  
 Pathways:Others  
 Receptor:Others  
 Boiling pt:549.1C  
 Melting pt:  
 Solubility:Soluble in DMSO  
 Appearance:White Powder  
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
 Stock solution can be stored below -20 °C for several month.

## References

1. Boissier JR, et al. Ann Pharm Fr. 1970 Jul-Aug;28(7):497-509.

Note: This product is for in vitro research use only