

## Hydrocortisone Valerate

Catalog No: #S0276


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Hydrocortisone Valerate
Brief Description	Inhibitors
Purification	98.00%
Target Name	GR inhibitor
Calculated MW	446.58
Formulation	C26H38O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:CCCCC(=O)O[C@@]1(CC[C@@H]2[C@@]1(CC[C@@H]([C@H]3[C@@H]2CCCC4=CC(=O)CC[C@]34C)O)C(=O)COPathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:595.1CMelting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Suzuki G, et al. Sci Total Environ. 2015 May 9;527-528C:328-334.

Note: This product is for in vitro research use only