

Orciprenaline sulfate

Catalog No: #S0277


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Orciprenaline sulfate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adrenergic Receptor antagonist
Calculated MW	260.3
Formulation	C ₁₁ H ₁₇ NO ₃ O ₅ H ₂ O ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory

systemSMILES:OS(O)(=O)=O.CC(C)NCC(O)C1=CC(O)=CC(O)=C1.CC(C)NCC(O)C1=CC(O)=CC(O)=C1Pathways:GPCR/G

ProteinReceptor:β-adrenergic receptorBoiling pt:417.5CMelting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Kimura M, et al. J Pharmacol Exp Ther. 1999 Oct;291(1):171-80.

Note: This product is for in vitro research use only