

## Escin

Catalog No: #S0278


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Escin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	1131.28
Formulation	C55H86O24
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism

systemSMILES:CC=C(C)C(=O)OC1C([C@@]2(C[C@@]3(C(=CCC4[C@]3(CCC5[C@@]4(CCC([C@]5(C)CO)OC6C(C(C(O6)C(=O)O)OC7C(C(C(O7)CO)O)O)O)OC8C(C(C(O8)CO)O)O)O)C)C)C2CC1(C)C)O)CO)OC(=O)CPathways:OthersReceptor:OthersBoiling pt:1140.6CMelting pt:Solubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Yu Z, Su P. Pharmazie. 2013 Jun;68(6):428-30.

Note: This product is for in vitro research use only