

## Metaraminol bitartrate

Catalog No: #S0279


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Metaraminol bitartrate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor agonist
Calculated MW	317.29
Formulation	C13H19NO8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:CancerSMILES:C[C@@H]([C@@H](C1=CC(=CC=C1)O)O)N.[C@@H]([C@H](C(=O)O)O)(C(=O)O)O

Pathways:GPCR/G  
 ProteinReceptor: $\alpha$ -adrenergic receptorBoiling pt:357.9CMelting pt:Solubility:Soluble in WaterAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Tatsuta M, et al. Int J Y. 1998 Jul 29;77(3):467-9.

Note: This product is for in vitro research use only