

## Sulfogaiacol

Catalog No: #S0284


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sulfogaiacol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	242.29
Formulation	C7H7KO5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:[K+].COC1=C(O)C=CC(=C1)S([O-])(=O)=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:34CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Armand C, et al. Therapie. 2004 Sep-Oct;59(5):547-53.

Note: This product is for in vitro research use only