

Bephenium Hydroxynaphthoate


Catalog No: #S0286

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bephenium Hydroxynaphthoate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	443.53
Formulation	C ₂₈ H ₂₉ NO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: OC(=O)C1=CC2=CC=CC=C2C=C1[O-].[C+](N)(C)CCOC1=CC=CC=C1CC1=CC=CC=C1 Pathways: Others Receptor: Others
 Boiling pt: Melting pt: 170-171°C Solubility: Soluble in DMSO Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. MEHROTRA MP, MALAVIYA US. Indian J Med Sci. 1963 Dec;17:930-3.

Note: This product is for in vitro research use only