

Procodazole


Catalog No: #S0288

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Procodazole
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	190.2
Formulation	C ₁₀ H ₁₀ N ₂ O ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer SMILES: OC(=O)CCC1=NC2=CC=CC=C2N1 Pathways: Others Receptor: Others Boiling pt: Melting pt: 229-231°C Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only