

Aminacrine

Catalog No: #S0289


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aminacrine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	194.23
Formulation	C13H10N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: NC1=C2C=CC=CC2=NC2=CC=CC=C12 Pathways: Others Receptor: Others Boiling pt: Melting pt: 238-240°C Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ju W, et al. *Retrovirology*. 2014 Jun 2;11:43.

Note: This product is for in vitro research use only