

Flurbiprofen

Catalog No: #S0291


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Flurbiprofen
Brief Description	Inhibitors
Purification	99.00%
Target Name	COX inhibitor; MRP inhibitor
Calculated MW	244.26
Formulation	C15H13FO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CC(C(O)=O)C1=CC(F)=C(C=C1)C1=CC=CC=C1
 Pathways: Immunology/Inflammation;
 Neuroscience
 Receptor: COX; MRP
 Boiling pt: Melting pt: 110-112°C
 Solubility: Ethanol: 49 mg/mL (200.6 mM); DMSO: 49 mg/mL (200.6 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Rieke CJ, et al. J Biol Chem. 1999 Jun 11; 274(24):17109-14.

Note: This product is for in vitro research use only