

Metyrapone


Catalog No: #S0293

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Metyrapone |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others inhibitor |
| Calculated MW | 226.27 |
| Formulation | C ₁₄ H ₁₄ N ₂ O |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area:Endocrine system SMILES:CC(C)(C(=O)C1=CN=CC=C1)C1=CN=CC=C1 Pathways:Others Receptor:Camphor 5-monooxygenase Boiling pt:384.4°C Melting pt:52-55°C Solubility:sparingly soluble in Water Appearance:White Crystal Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Shiro Y, et al. Biochemistry. 1995 Jul 18;34(28):9052-8.

Note: This product is for in vitro research use only