

## Digitoxin

Catalog No: #S0295


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Digitoxin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	764.94
Formulation	C41H64O13
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular

SystemSMILES:C[C@@H]1[C@H]([C@H](C[C@@H](O1)O[C@@H]2[C@H](O[C@H](C[C@@H]2O)O[C@@H]3[C@H](O[C@H](C[C@@H]3O)O[C@@H]4CC[C@]5([C@@H](C4)CC[C@@H]6[C@@H]5CC[C@]7([C@@]6(CC[C@@H]7C8=CC(=O)OC8)O)C)C)O)OPathways:OthersReceptor:organic anion transporterBoiling pt:902.3CMelting pt:240CSolubility:Water: 3.9 mg/lAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Bossuyt X, et al. J Pharmacol Exp Ther. 1996 Mar;276(3):891-6.

Note: This product is for in vitro research use only